



ZAPR™ Red Blood Cell Lysing Buffer

Product Name: ZAPR™
Product Use: Red Blood Cell (RBC) Lysing Buffer solution used for tissue and cell processing
Product Codes: ZAPR-100
Features: GMP quality; sterile; USP grade materials; no animal components; no human proteins; no antibiotics



Manufacturer: INCELL Corporation LLC
Address: 12734 Cimarron Path
City/State/Zip: San Antonio, TX 78249
Phone: 800.364.1765; 210.877.0100
Email: info@incell.com
FAX: 210.877.0200

Distributor: INCELL Corporation LLC
Address: 12734 Cimarron Path
City/State/Zip: San Antonio, TX 78249
Phone: 800.364.1765; 210.877.0100
Email: info@incell.com
FAX: 210.877.0200

General Description

ZAPR™ is INCELL's Red Blood Cell (RBC) Lysing Buffer Solution. ZAPR™ is sterile, "xeno-free", GMP manufactured, chemically defined and ready-to-use. All components used to manufacture ZAPR™ are USP grade. ZAPR™ has a Master File with the FDA. Red blood cells (RBCs) are commonly found in processed tissues since they are a blood component. ZAPR™ is used as a rinsing solution to clear RBCs from isolated tissue samples and from products intended for clinical use. ZAPR™ was developed for use with human tissues, but it is also an effective RBC lysing buffer that is effective with other mammalian source tissues.

Formulation and Packaging

GMP manufactured ZAPR™ RBC lysis buffer is formulated with USP grade components. ZAPR™ is available in 100 mL square PETG media bottles. Product or catalog number is ZAPR-100. Larger volumes or custom packaging (e.g., bags) can be special-ordered with 3 to 4 weeks lead time and purchase of at least 100 units.

Use and Methods

ZAPR™ is used as a post processing rinse to lyse residual RBC in processed tissue or cell samples, including cells from blood. Pre-rinse the collected tissues to minimize blood contamination. Process tissues by tissue-appropriate methods (e.g., mincing, enzymatic digestion, chelating agents, rinsing) to release and separate cells. Centrifuge cell suspension to pellet the cells at a g force appropriate for the cells being released. A density gradient buffer may be used to fractionate some cell suspensions, if appropriate, particularly for blood or bloody samples. After the suspension processing solution or buffer is removed, add 5 volumes of ZAPR™ a cell pellet or fractions containing a smaller percentage (<20%) of RBCs in the cell population. For whole blood, use 10 volumes of ZAPR™ per volume whole blood (e.g., use 100 mL ZAPR™ per 10 mL whole blood).

The ZAPR™-cell suspension is gently mixed to ensure homogeneity and to minimize RBC clumping. The suspension is left at room temp (15°C to 30°C) for 10 minutes. The suspension is diluted at least 2-fold with cell appropriate media or buffer, then again centrifuged to gently pellet the cells as per the protocol for that cell type. If the pellet is still red due to RBCs then repeat the ZAPR™ treatment steps again and do the cell pipetting gently because the cells may be somewhat fragile. Re-suspend the cells in a storage media, storage solution or buffer, serum, or culture media depending on intended use.

Manufacturing

ZAPR™ Red Blood Cell Lysing Buffer is manufactured by sterile 0.22 µm filtration and packaging, using cGMP standards in an ISO Class 7 clean room and ISO Class 5 biosafety cabinet. Raw materials are pre-tested, and the final product is assayed for endotoxin, sterility (bacteria, fungi, and mycoplasma) and other quality specifications and criteria (see Table 1) prior to release and preparation of the Certificate of Analysis (COA) for the lot.

ZAPR™
RBC Lysing Buffer



Table 1. Data Collected for the COA
Specifications **Acceptance Criteria**

| | |
|------------|-----------------------------|
| pH | 7.2 to 7.6 |
| Osmolality | 260 to 320 mmol/kg |
| Endotoxin | <0.1 EU/mL |
| Mycoplasma | None detected |
| Sterility | No microbial growth |
| Visual | Colorless, clear solution |
| Expiration | 24 months after manufacture |

Storage of ZAPR™ RBC Lysing Buffer

ZAPR™ Red Blood Cell Lysing Buffer solution can be stored at room temperature (15° to 30°C). The shelf-life is 24 months from the date of manufacturing.

Table 1. Examples of Human and Animal Tissues, Cells and Other Biological Products for Potential Use of ZAPR™ Suitable for Bio-Processing

| Tissue Types | Primary Cell Type(s) Derived |
|--|--|
| Adipose (Fat) | Mesenchymal Stem Cells; Stromal vascular fraction regenerative cells; adipose cells |
| Bone Marrow Bone | Hematopoietic and mesenchymal stem cells; various types of renewable progenitor cells; Endothelial cells; entire population; bone & cartilage progenitors; other cells |
| Brain & Neural (Spinal) | Progenitors; Induction of differentiation |
| Breast | Primary epithelial cells and/or mesenchymal support cells |
| Cartilage | Chondrocytes & Progenitors |
| Colon; Gastrointestinal; oral; esophageal | Primary epithelial and/or mesenchymal support cells or complex tissues in organ-like cultures |
| Dental Pulp | Hematopoietic and mesenchymal stem cells; various types of progenitor cells; similar to bone marrow |
| Heart | Primary muscle and mesenchymal cells |
| Kidney | Primary epithelial cells and/or mesenchymal support cells |
| Liver | Primary epithelial cells and/or mesenchymal support cells |
| Lung | Primary epithelial cells and/or mesenchymal support cells |
| Lymph Node | Lymphocytes and mesenchymal support cells |
| Muscles (Peripheral; Heart; Smooth) | Pericytes; Progenitor muscle cells; Mesenchymal or Stromal Stem Cells; other regenerative cells |
| Pancreas; Other Neuroendocrine organs | Pancreatic islet beta and acinar cells; other organs (e.g., adrenal) |
| Parathyroid | Adenoma cells isolated and stored for potential autologous re-transplantation |
| Peripheral or apheresis blood | Various types of white blood cells, including lymphocytes, macrophages, etc.; Circulating or mesenchymal cells; endothelial cells. |
| Placenta and birth tissues | Trophoblasts; Syncytiotrophoblasts; Endothelial cells; Hematopoietic and mesenchymal stem cells; various renewable progenitor cells from amnion and chorion |
| Prostate | Primary epithelial cells and/or mesenchymal support cells |
| Skin (adult; newborn) | Epidermal keratinocytes; biopsies; hair follicles Dermal Fibroblasts; Mesenchymal cells |
| Spleen | Lymphocytes and mesenchymal support cells |
| Tumors | Epithelial, mesenchymal, lymphoid |
| Umbilical cord | Hematopoietic and mesenchymal stem cells; various types of renewable progenitor cells; Endothelial cells |

*Use of ZAPR™ RBC Lysis solution is based on protocol and processing needs.

Master Files Applications Note

INCELL media and solutions are in FDA Drug and Device Master Files but have not been tested by INCELL for any specific diagnostic or therapeutic use. To request use of a Master File call, FAX, or email to info@incell.com.

Animal Component Free. INCELL certifies that the Product is “animal component free” per the following criteria:

- No animal derived supplements are added by INCELL or come from raw materials supplied as components of the Product.
- Product does not come into contact with animal derived material during manufacturing, processing, handling, or packaging.
- This certification applies only to the condition of the above-described Product in its unopened package, and INCELL assumes no responsibility for a Product failing to meet this Statement after handling or use after opening the package.
- Signed “Animal Component Free Statement” can be provided on request.

Ordering: Contact INCELL Corporation

Toll Free: 800.364.1765
info@incell.com
FAX: 210.877.0200

Technical Assistance

The scientists at INCELL are available to discuss the media or special needs of your cells, and to assist you in your cell culture applications. Call 1-800-364-1765 or e-mail info@incell.com.